


# Medcold PG300 Pharmacy Fridge

- Temp range: +2°C/+8°C
- External digital controller accurate to 0.1°C
- Forced air cooling
- Electronic controller
- Max & min temperature recording
- High & low alarm
- Visual & Audible alarm
- Glass door
- Lockable



PG300

## Special Features

- Safety thermostat to prevent freezing of medicaments
- Door alarm
- LED low energy internal lighting
- Auxiliary output for remote alarm / BMS system

Model	Dimensions ( L x D x H mm )	Capacity ( L )	Shelves
PG300	595 x 608 x 1563	278	6
Cool Tray Drawer*	457 x 404 x 60	-	-

\*with 3 lanes & 12 internal divisions - Out of 2 years warranty

Please Refer to Manufacturers Catalogue for Full Product Specifications