

Medline MED 180G Pharmacy Fridge

- Temp range: +2°C/+8°C
- External digital controller accurate to 0.1°C
- Forced air cooling
- Electronic controller
- Max & min temperature recording
- High & low alarm
- Visual & Audible alarm
- Glass door
- Lockable
- LED low energy internal lighting
- Auxilary output for remote alarm/ BMS system
- Door alarm
- Warm Air Dispersal System



- Designed to be positioned undercounter with limited ventilation/ air flow.
- Warm Air Dispersal System, extracts excess heat from the rear.
- Increased fridge lifespan due to warm air dispersal system.

MED 180G

Model	Dimensions (L x D x H mm)	Shelves	Capacity (L)	Weight (Kg)	Price
MED 180	600 x 600 x 850	3	150		

Please refer to manufacturers catalogue book for full product specifications and exact product colours